Thopaz™

Thoracic Drainage System

- EN Instructions for use
- DE Gebrauchsanweisung
- FR Mode d'emploi
- IT Istruzioni per l'uso
- NL Gebruiksaanwijzing
- SV Bruksanvisning
- DA Brugsanvisning
- ES Instrucciones de uso
- PT Instruções de utilização
- PL Instrukcja obsługi

-10 kPa / 5 L/min

Technical specifications – Technische Daten – Données techniques –
Specifiche tecniche – Technische specificaties – Tehniska specifikatio-
ner – Tekniske specifikationer – Especificaciones técnicas – Especificações
técnicas – Specyfikacja techniczna

Low vacuum

- 10 kPa
- 75 mmHg
- 100 mbar/cm H₂O

H x W x D

95 x 170 x 235 mm

3.74 x 6.69 x 9.25 inch

low flow

5 L/min.

- 42.5 dB(A) 1 L -2.5 kPa

Max. noise level
Max. Lautstärke
Niveau de bruit max.
Livello di rumore max.
Max. geluidsniveau

Max. ljudnivå
Max. støyniveau
Nível de ruído máx.
Nível de ruído máximo
Maks. poziom głośności

1 kg
2,2 lbs

Without canister
Ohne Behälter
Sans récipient
Senza contenitore
Zonder opvangpot

Utan behållare
Uden beholder
Sin recipiente
Sem recipiente
Bez pojemnika

78 dB(A)

Alarm noise level
Alarm Lautstärke
Alarme de niveau de bruit
Livello di rumore allarme
Geluidniveau alarm

Larm ljudnivå
Alarmernes støyniveau
Nível de ruído de la alarma
Nível de ruído de ruído
Poziom głośności alarmu

IP33

ISO 9001
ISO 13485
CE (93/42/EEC), IIa

CE
0123

Operation
Betrieb
Fonctionnement
Funzionamento
Gebruik

Drift
Drift
Foncionamiento
Funzionamento
Działanie

Operation
Betrieb
Fonctionnement
Funzionamento
Gebruik

Drift
Drift
Foncionamiento
Funzionamento
Działanie

Operation
Betrieb
Fonctionnement
Funzionamento
Gebruik

Drift
Drift
Foncionamiento
Funzionamento
Działanie

Transport/Storage
Transport/Lagerung
Transport/Stockage
Trasporto/Conservazione
Transport/Opslag

Transport/Förvaring
Transport/Opbevaring
Transporte/Alcanamiento
Transporte/Amazamento
Transport/Magazynowanie

Transport/Storage
Transport/Lagerung
Transport/Stockage
Trasporto/Conservazione
Transport/Opslag

Transport/Förvaring
Transport/Opbevaring
Transporte/Alcanamiento
Transporte/Amazamento
Transport/Magazynowanie

20 W
12 VDC

Li-Ion

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Caution: Federal US law restricts this device to sale by or on the order of a health professional.

Subject to change - Änderungen vorbehalten - Sous réserve de modifications - Soggetto a modifica - Wijzigingen voorbehouden - Med reservation för förändringar - Med forbehold for ændringer - Sujeto a cambios - Sujeito a alteração - Može ulec zmianie

I Congratulations

With Thopaz, you have acquired an innovative digital thoracic drainage system that sets new standards in therapy and thoracic drainage management.

Thopaz has an electronic measuring and monitoring system with optical and acoustic status indications. The device is a dry system, which means no fluids are necessary for operation. Important information regarding the course of therapy is indicated digitally and graphically on the display. Upon completion of therapy, it can be transferred to a PC.

The handy size and extreme quietness of operation are very pleasant for the patient as well as those in the vicinity.

II Intended use

Thopaz can be used for the following purposes:

- to remove air and/or fluids which gather in the intrapleural space as a result of illness, injury or surgery;
- to regain or maintain negative pressure in the intrapleural space and to unfold the lungs.

Thopaz is intended for use on patients in all stages of the hospital. The unit should only be operated by properly trained staff. These persons must not be hard of hearing or deaf and must have adequate visual faculty.

Thopaz is portable and, thanks to a rechargeable battery, can be operated independently of the mains supply. The built-in electronics allow you to make the necessary adjustments and to track the effectiveness of the therapy during operation. Acoustic and optical signals are triggered by variances during the course of treatment.

Indications

The system is used when negative pressure is required in combination with a thoracic drainage catheter.

Contraindications

Thopaz cannot be used if:

- a pressure greater than -10 kPa or a flow of more than 5 L/min is required;
- applying active suction would put the patient at risk.

Medela does not recommend connecting bilateral thoracic drains to one Thopaz unit. In such cases, the use of two Thopaz units is recommended.

III Warnings / Safety instructions

Thopaz is approved exclusively for use as described in these instructions. Medela can only guarantee the safe functioning of the system when Thopaz is used in combination with the original Medela accessories (canisters, tubing sets, carrying strap, mains adapter, docking station).

Thopaz is EMC-tested in conformity with the requirements of IEC 60601-1-2:2007 and can be used in the vicinity of other EMC-tested devices that fulfil the requirements as outlined in the IEC 60661-1-2 standard. Untested HF sources, radio networks or the like can impair the function of the device and should not be operated in the vicinity of Thopaz.

Please read and observe these warnings and safety instructions before putting into operation.

Warnings

- Before connecting Thopaz to the power supply, please verify that the supply voltage corresponds to that given on the device specification plate.
- Data transfer via USB is not possible while the pump is connected to the patient.
- The device is not suitable for use while bathing, showering or in a hazardous explosive environment.
- Do not use Thopaz in MRT (Magnetic Resonance Tomography).
- Do not dry Thopaz with microwaves.

Safety instructions

- Thopaz must remain in an upright position during use.
- The patient should be regularly monitored according to internal hospital guidelines.
- The pressure range to be set must be determined by a physician in accordance with the age and weight of the patient.
- The secretions gathered must be monitored more regularly or frequently in accordance with the age and weight of the patient.
- Thopaz is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect Thopaz.
- A replacement device must always be available for patients for whom a breakdown of the device can lead to a critical situation.
- In none of the following cases may Thopaz be placed in operation:
 - if the power cord or plug is damaged
 - if the device is not functioning properly
 - if the device is damaged
 - if the device has apparent safety defects.
- Keep the power cord away from hot surfaces.
- Keep the AC adapter connector away from moisture.
- Disconnection from the power supply is achieved by unplugging the mains adapter.


These instructions for use must be kept for later reference.

IV Initial startup

Check the condition of the delivery package and that it is complete.

Thopaz™
079.0000
079.0002
079.0003

Mains adapter
077.0148

Thopaz
Instructions for use
200.0685
200.0686
200.2004


Thopaz CD
079.0035

Do not switch Thopaz on before first charging the battery!


Initial startup

1. Select plug

2. Charging battery for the first time


Connect Thopaz to mains

► *The battery must be fully charged before first usage.*

- a) **Disconnect Thopaz from mains power** – pull on the plug housing.

Do not pull on the cable or the bend protection!

- b) Thopaz switches off automatically.




3. Select language

- a) Switch Thopaz on with [].
- b) Select language (follow instructions on display), confirm with 'OK'.
- c) Switch Thopaz off with []. Thopaz is ready for use.

► *If Thopaz was switched on before the battery was fully charged (step 2), follow the instructions on the display.*

► *Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended to acquire the reference values.*


V Overview of Thopaz

VI Disposables


Canister with/without solidifier 0.3 L and 0.8 L

Material: Polypropylene
 Reading accuracy: +/- 2.5% (in upright position)
 Solidifier: 0.3 L = 12 g / 0.8 L = 28 g

Tubing sets single / double

Material: PVC (medical grade)
 sterile, twofold packaging
 Ø / Length: 5 mm / 1.5 m

► *Overflow protection / bacteria filter see 'Additional information'.*

VII Preparation for use

Use only after instruction by qualified personnel.


► 'Pressure' generally implies 'negative pressure'.

1

Check required parts

- Thopaz
- Tubing set (single / double)
- Canister
(usable capacity in upright position
0.3 L / 0.8 L)

without solidifier

with solidifier


2


Connect tubing set

Open external packaging.

Insert the smaller of the two connectors horizontally in the direction of the arrow.

3


Snap in canister

Unpack canister, position the openings at the top and attach the bottom part of the canister to the pump.

► *Push the canister into Thopaz until you hear a click.*

4

Switching on

1. **Switch on Thopaz** with [] – self-test starts (check acoustic signal – [beep]).


 *Do not switch Thopaz on if the patient is already connected.*

If the self-test is not successful, see the troubleshooting instructions on the display (Chapter XII).

2. **New patient? yes/no**

- **'yes'** means that Thopaz will issue a new therapy number (recommended for a new patient).
- **'no'** means that the therapy number remains unchanged (recommended for the continued treatment of the same patient). The graphic presentation is retained.

Therapy number

 *The therapy number is important for data transfer to the PC.*


5

Carry out functional check

1. **Seal off the conical connector of the patient tubing with your thumb** (through the packaging).

 *The tubing must not be kinked as otherwise the measuring tube will be closed.*

2. • **Switch on pressure** by pressing 'on' – pressure is built up.

- **Check flow value**

 *Flow decreases:*

Thopaz is ready for use

Flow does not decrease:

*Leak in the system;
check tubing connections,
replace canister and/or
tubing*

3. **Switch to Standby mode:** Press 'Standby' (press for > 3 sec.)


4. **Thopaz is in Standby mode and is ready for use with the factory settings.**

1

Connect patient catheter in accordance with the physician's directions.

► Remove the inner packaging from the tubing.

2

Usage

1. **Switch on Thopaz** with [] – self-test starts.

► Do not switch Thopaz on if the patient is already connected.

2. **New patient? yes/no.**

3. • **Switch on pressure** by pressing 'on' – the pre-selected pressure is built up and kept constant.

- **Check flow value.**

► The acoustic alarm function "Leak in system - 301" is deactivated for approximately 5 minutes after switching on.


3

Check therapy progress

There are two different modes available:

Data mode and **Graph mode**


In the data mode, the set pressure value and the current flow value (equivalent to parenchymal leakage) are digitally displayed.

Data mode

Flow values between 0 and 1000 mL/min = **10 mL steps**
Flow values over 1000 mL/min = **100 mL steps**

► Lowest displayed flow value in data mode = 10 mL/min.


Graph mode

With 'Graph', the Graph mode is activated and the flow and pressure progression are displayed on a time line. After 60 seconds, the display switches back to the data mode.

With 'Data', you can switch to the Data mode at any time.

Graphic presentation

The Graph mode shows the flow and pressure progression as a function of time. The graph is updated every ten minutes. The displayed graph is first meaningful after 4 hours.



- The current flow and pressure values are indicated in the information field.
- The left scale shows the pressure in the desired unit (max. -10 kPa). The pressure is presented as a filled area.
- The right scale shows the flow in the desired unit (max. 1000 mL). The flow values are presented as a line.

Zoom mode and catheter check

1. The zoom function is activated with the **selection button** []. The flow scale is reduced to 100 mL/min.

Return to the Graph mode with the **selection button** [].

2. In zoom mode, the pressure difference between inspiration and expiration is displayed at a flow of 0 mL/min. This pressure difference is an indication for the respiration cycle of the patient and confirms that the catheter is not clogged.

catheter open



catheter clogged?

Indicator only visible in zoom mode at a flow of 0 mL/min

If Thopaz does not show the respiratory movements, request the patient to breathe slowly and deeply.




Switching off

1. Close patient tubing with the tube clamp.
2. Switch off pressure by pressing 'Standby' (press for > 3 sec.)
– Thopaz switches to standby mode.
3. Release, remove and seal the canister.
4. Dispose of canister and patient tubing in accordance with internal hospital guidelines.
5. Switch off Thopaz [].

Option: Transfer data to PC – see Chapter XI.

IX Replace canister

The canister is replaced on the basis of a visual check or according to the instructions on the display (warning signal).


1. Prepare sterile canister (with/without solidifier).
2. Close patient tubing with the tube clamp.
3. **'Standby'** (press for > 3 sec.)
Thopaz switches to the Standby mode.
4. Release and remove canister.
5. Unpack new canister, position into Thopaz and snap in (Chapter VII).
6. Press **'on'**. Pressure is built up.
7. Open the tube clamp.
8. Check flow value.
9. Seal the used canister with the sealing caps.

Only for canisters with solidifier:

a) Press – to open the solidifier chamber.

b) Shake – to activate the solidification process.


10. Dispose of the used canister in accordance with internal hospital guidelines.

X Change settings

▶ The settings may only be changed by the physician or on medical orders.

Operating mode

Change pressure

- Press [▶] and [◀] simultaneously. The pressure is indicated on the display.
- With the **selection buttons** set the desired pressure and confirm with 'OK'.
- To activate the **gravity mode** press the [▶] button and confirm with 'OK'.

▶ For patients who are to be treated by gravity drainage (= water seal), the gravity mode can be activated. This mode corresponds to a pressure of -0.8 kPa / -8 cmH₂O / -8 mbar / -6 mmHg.

Standby mode

Making settings

- With the **selection buttons**, select the desired parameter and confirm with 'OK'.
- With the **selection buttons**, select the desired setting and confirm with 'OK'.
- With the 'back' button, the settings mode is ended and the standby mode is activated.


Settings 1/2		Selection
Pressure	- 2.0	max. -10 kPa / -100 cmH ₂ O
Unit Pressure	kPa	kPa / cmH ₂ O / mbar / mmHg
Unit Flow	mL/min	mL/min / L/min
Leakage	off	off / on
Set as standard	no	no / yes
back	OK	

► *cmH₂O is not a SI unit according to EU Directive 80/181/EEC.*

Settings 2/2		Selection
Language	English	Languages according to list
Version	1.2x	Display
Thopaz Number	1234567	Display
Runtime	175 h	Display
back	OK	

Set as standard

Thopaz is delivered with the following factory settings.

These settings can be changed and saved as the new standard and are active whenever Thopaz is switched on.

To save the selected settings, confirm parameter 'Set as standard' with 'yes'. Press 'back' to activate the standby mode.

XI Data transfer to PC


With the software ThopEasy all stored data can be transferred to a PC as documentation and for completion of the patient records. The data can be expanded with patient information, saved and printed.

The operating instructions for data transfer and the software ThopEasy can be found on the Thopaz CD included in the scope of delivery.

XII Alarms

Thopaz distinguishes between warnings, alarms and errors.

If Thopaz detects any of these situations, an acoustic warning signal sounds and a description of the problem appears on the display. By pressing the two selection buttons simultaneously, the acoustic alarm is suppressed for 60 seconds and the troubleshooting instructions appear on the display.

	Visible on display				
	Error number	Problem description	Troubleshooting instructions	Remarks / Potential source of error	Pressure
Warning	401	Battery low	Charge battery – connect Thopaz to mains	Remaining time of battery approx. 30 min.	✓
	306	Canister full	Replace canister	If the canister is not full, see Chapter XIII	✓
	405	Standby mode	Switch Thopaz on or off	In standby mode after 5 minutes	✗
	404	Leakage	• Check patient status • Check system for leaks • Continue with 'OK'	• Irregularity in the course of therapy • Indication of leakage • See Chapter XIII	✓
	402	USB connection not permitted	Unplug USB cable	USB connection not permitted during operation or in the vicinity of the patient	✓
	406	Temperature too high	Prepare backup device	Cool Thopaz (for example with a fan, hair dryer, etc.). Do not expose Thopaz to direct sunlight.	✓
Alarm	305	Battery empty	Charge battery – connect Thopaz to mains		✗
	302	System clogged	• Check that tubing is clear/not twisted • Replace canister	• Tubing is kinked or clogged • Filter in canister clogged	✗
	301	Leak in system	Check system for leaks	• Disconnection • Missing seal between Thopaz and canister	✗
	311	Selftest failed	Snap the canister out and in again	Do not switch Thopaz on if the patient is already connected.	✗
	313	Filter clogged	• Replace canister • Continue with 'on'	Filter in canister clogged	✗
	315	Thopaz overheated	Replace Thopaz		✗
Internal error		Internal error XXX	1 Switch Thopaz off and on again 2 Inform Medela Customer Service	If error message repeated, inform Medela Customer Service	

XIII Additional information

Regulation

The pressure is regularly checked and adjusted by Thopaz. The set pressure corresponds to the pressure of the entire system including the patient.

Safety pressure limit

The adjustable pressure range for Thopaz is 0 to -10 kPa.

Pressure value	cmH ₂ O	mbar	mmHg
-7 kPa	-70	-70	-53
-10 kPa	-100	-100	-75

If pressure values greater than -7 kPa are set, the following warning appears in the display:
Too high pressure can cause pain and injuries to the patient.

This warning must be acknowledged with '**OK**' before the pressure can be increased further.

Tube flushing

Flushing of the tubing occurs every 5 minutes or when Thopaz detects a siphon and prevents clogging of the patient tubing.

Safety chamber in canister

For appropriate, safe operation Thopaz must be in an upright position. If Thopaz tips, the special construction of the safety chamber in the upper region of the canister protects the hydrophobic filter from premature clogging. The pressure is maintained.

Important: immediately set Thopaz upright again.

This function is disabled if the maximum filling capacity of the canister is exceeded.

Canister full

If secretions line the walls of the canister, the warning 'canister full' may be activated prematurely. The warning 'canister full' is triggered only once. After cancelling the warning with the selection buttons, it is deactivated and a symbol [] flashes in the run mode.

After deactivation of the warning, the user is responsible for checking the canister filling level. If the hydrophobic overflow/bacteria filter gets clogged by fluids, the alarm 'System clogged' is triggered.

The warning is reactivated under the following conditions:

- if Thopaz is switched off and then on again – for example, on replacing the canister.
- if the level sensor is not activated for at least 10 minutes.

Overflow protection / bacteria filter

A hydrophobic overflow protection / bacteria filter in the canister as well as in the tube sets (in the connector to the pump) protects Thopaz from the entry of fluids and protects the environment from contamination.


Leakage

A significant and sustained increase in the flow rate indicates that there is a leakage in the system or on irregularity in the course of therapy. This function can be switched on as required - see Chapter X.

XIV Optional accessories

079.0031
Carrying strap

079.0037
Docking station
(is delivered without mains
adapter)

079.0034
USB - standard cable
(not included in scope of delivery)

XV Cleaning

Thopaz with accessories should be cleaned / disinfected after every use.

General notes

- The internal hospital guidelines take precedence.
- Disconnect Thopaz from the mains before cleaning / disinfection.
- Classification according to RKI (Robert Koch Institut) guidelines: non-critical.
- Wear suitable gloves for cleaning / disinfection.
- Dispose of fluids such as blood and secretions and the parts contaminated with them in accordance with internal hospital guidelines.
- Detailed information according to Medela cleaning instructions Art. No.200.2391.

	Cleaning	Disinfection	Rinse cycle in washing machines	Sterilisation
Thopaz / Docking station Accessories <i>(e.g. docking station, carrying strap)</i>	 Wipe off with a damp cloth.	 Wipe off with a disinfecting agent, e.g., Mikrozid AF Liquid (supplier Schülke + Mayr).	 • Not sterilisable. • Rinsing in a washing machine is not permitted. • Immersion disinfection, thermal disinfection and ultrasound cleaning are not permitted.	
Disposables <i>(canister, tubing set)</i>		Caution: This is a single use product not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and / or biological characteristics. Reuse could cause cross contamination.		

Thopaz can be disinfected with the disinfecting agent group 'alcohol'.
Do not use other cleaning agents (e.g. Terralin) as they could damage the plastic housing.

XVI Power supply / Battery management

Thopaz can be operated from the mains power supply or with the integrated lithium-ion battery. The battery is charged during mains operation.

The battery duration is dependent on the run-time of the pump. This is influenced by the extent of parenchymal leakage and the set pressure. Thopaz does not run continuously but only switches on when the actual and nominal values differ.

During continuous operation, Medela guarantees a minimum run-time of 4 hours in battery operation. This case will not occur in practice, which results in real battery run-times of > 10 hours.

If Thopaz is not used very often, the battery must be charged approximately every 6 months to ensure optimum functioning.

Backup battery

If an internal defect occurs (broken cable, battery defective), an acoustic warning signal is sounded for at least 3 minutes. Power is supplied by the backup battery. Replace Thopaz immediately.

XVII Guarantee

Thopaz is guaranteed for 2 years after the date of delivery when used in accordance with these instructions. The manufacturer is not liable for any damage or consequential damage caused by incorrect operation, inappropriate use as well as use by unauthorised persons.

XVIII Routine check / Safety test

Each time the device is switched on, a self-test is carried out which checks the internal functions of Thopaz. A "beep" indicates that Thopaz is switched on. Safe operation is ensured by the functional check including canister and tubing set to be carried out before each use. No further maintenance or routine checks are necessary. The safety instructions must be observed.

Medela recommends subjecting Thopaz to a safety check once a year. The work to be carried out is listed in Appendix B in these instructions.

XIX Disposal

Thopaz comprises metals and plastics and should be disposed of in accordance with the European directives 2002/95/EC and 2002/96/EC. Additional, local guidelines must also be observed.

079.0000 Thopaz drainage pump with Instructions for use 200.0685

079.0002 Thopaz drainage pump with Instructions for use 200.0686

079.0003 Thopaz drainage pump with Instructions for use 200.2004

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Zmiana kanistra 202, 206

Electromagnetic emissions

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such an environment.

Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 60000-3-3	Complies	

Warning – Thopaz should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz should be observed in order to verify normal operation in the configuration in which it will be used.

Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that they are used in such an environment.

Immunity tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	< 5% UT (> 95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz requires continued operation during power mains interruptions, it is recommended that Thopaz be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) Magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

Electromagnetic immunity

Thopaz is intended for use in the electromagnetic environment specified below. The customer or the user of Thopaz should assure that Thopaz is used in such an environment.

Immunity Tests	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of Thopaz, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz is used exceeds the applicable RF compliance level above, Thopaz should be observed in order to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distance between portable and mobile RF communications equipment and Thopaz			
Thopaz is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Thopaz can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and Thopaz, as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter	Separation distance according to frequency of transmitter		
W	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

B Safety test

EN IEC 62353: 2008-8 Recurrent Test and Test after Repair

For safety tests, it is assumed that Thopaz is maintained and repaired throughout its service life strictly and exclusively by Medela authorised service centres.

For Thopaz as a device in protection class II (EN IEC 60601-1), the safety tests are confined to visual inspection of the housing and mains power cord for damage. This check must be performed prior to each use.

Devices of protection class II do not have a protective earth conductor; there is therefore no need to check the earth leakage current!

The Thopaz enclosure is made entirely of insulating material. Tests of the enclosure leakage current using common measuring instruments will therefore not reveal measurable values.

The applied parts of Thopaz are connected to Thopaz via non-conductive vacuum hoses, collection jars and intermediate tubings.


Even when suctioning a conductive fluid until the overflow protection device activates (9g/L NaCl; EN ISO 10079-1;1999), measurements of the patient leakage current using common measuring instruments will not reveal measurable values.

Thopaz pumps do not have patient circuits and functional earth connections.

Instructions for safety test:

▶ Patient must not be connected!
Thopaz must be disconnected from the mains (battery operation).

1. Start safety test (press buttons simultaneously).

2. Follow the instructions on the display.

Additional tubing
(e.g. $\varnothing$ 12x7 mm)
(not included in the
scope of delivery)

B Thopaz test form

DIN EN 62353: 2008-8

Hospital	Thopaz serial number:
	Firmware version:
Location:	Run time (Settings 2/2): h

Visual checks	NO	YES
Cracks, brittle areas, discolouration on housing and display		
Covers for ports are present and can be clicked in		
Standard rail and bed holders are present and can be opened		
Seal is present, is clean and not brittle		
Release button can be pressed, canister can be clicked in		
Thopaz is clean with no dirty areas		

Tests	OK	Not OK
Safety test Carry out test as described in Instructions for use. Thopaz must be disconnected from the mains (battery operation). Required test equipment: Vacuum gauge ~0-10 kPa Class 2.5% or better (E.G.077.1456)	☐	☐
Acoustic signal (beep) Audible?	☐	☐
Battery Connect Thopaz to mains adapter or place it in docking station.	☐ Battery is being charged (refer to Chapter IV)	☐ Battery is not being charged

Overall evaluation	Errors	Actions
No safety or functional technical errors were detected.	☐	Thopaz is ready for use
Errors that are not safety relevant were detected. No direct risk.	☐ 1 _____ 2 _____ 3 _____	Thopaz is ready for use, errors must be corrected by ____ / ____ / ____
Safety relevant errors were detected.	☐ 1 _____ 2 _____ 3 _____	Do not use Thopaz until the errors are corrected. Send Thopaz to Medela Service Centre.

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Medela AG, Medical Technology
Medela UK Ltd.
Medela Inc.